



ATTENTION: HANDLING INSTRUCTIONS

Please follow these instructions and handle this chair with care. Its owner's ability to live their life depends on its safe arrival to their destination. Thank you!

- **NOTE THE WEIGHT OF THIS CHAIR BELOW.**
It will require at least two people to lift it.
- **DO NOT** lift by the arm rests or any piece other than the base of the chair. Other pieces will break if used for lifting.
- **DO NOT** store on its side. Keep the chair upright.
- **ENSURE THE CHAIR IS IN "MANUAL" MODE WHEN PUSHING IT.**
 - ▶ **TO PUT THE CHAIR IN "MANUAL" MODE**, turn off the chair, and unlock each wheel. "Lock" the wheels when parked.
 - ▶ **DO NOT DRAG THE CHAIR.** If the wheels aren't turning easily when you push/pull it, this will result in damage to the chair.
- **FURTHER INSTRUCTIONS:**

INFORMATION AND CONTACT:

Owner:

Phone:

Weight:

Battery Type:

ATTENTION: HANDLING INSTRUCTIONS